



CareGard Quick Reference

Press the Volume button repeatedly until the volume level is comfortable.

Power/Home/Away lights up when AC Power is on.

Press it to turn Activity Check and Medication Reminders off when you are Away and back on when you are Home. The button flashes while you are Away.

Press the CALL button or your personal help button to signal for help if you having an emergency.

If the Reset/Status button is flashing, press it and the panel will tell you the current status. If there is a problem, consult your User Guide.

For more detailed information on the operation of the CareGard system refer to the CareGard User Guide.

